


EPFL-TTO EPFL Innovation Park J CH-1015 Lausanne Switzerland +41 21 693 70 23 http://tto.epfl.ch/ Research Contact

Carlotta Guiducci +41 21 693 78 13 carlotta.guiducci@epfl.ch TTO Contact Natalia Giovannini +41 21 693 35 90 natalia.giovannini@epfl.ch

Licensing Opportunity

TTO - Technology Transfer Office

Sharp edged electrodes for DLD devices



Ref. Nr 6.2106

Kevwords

Particle sorting, sample purification, microfluidics, dielectrophoretic particle manipulation

Intellectual Property

PCT/IB2020/059823 Priority 19.10.2020

Publications

DOI: 10.1002/admt.20190033

Date

28/01/2021

Sharp edged electrodes help to reduce the sorting size of DEP-DLD devices by creating high field regions where different sized particle streams split up

Description

Device, system and method for sorting particles comprised into a fluid that improve the current settings of deterministic lateral displacement (DLD)based microfluidic devices, so to sort small particles such as extracellular vesicles, while exploiting all the advantages and of well-known DLD throughput and dielectrophoresis techniques.

The device is characterized by the presence of plurality of three-dimensional а electrodes having a longitudinal section with a narrower first end and a wider, opposed second end. This creates high electric field gradients, where particles of different sizes follow different trajectories. gradient The high field enables an effective tuning of the sorting and enables to reach size nano size (currently down to 100 nm of size) sorting of particles suspended in liquids.

The project is currently supported by the enable grant of the Technology Transfer Office. We are currently working on the proof of concept for exosome purification.

Advantages

Combining mechanical and electric sorting $\ _{\bullet}$ technology enables a

- Sorting base on dielectric properties
- Increase throughput

The sharp edged electrodes enable

- More efficient sorting compared to round shaped electrodes, allowing to operate the device less than 20 V reaching 100 nm sorting size
- Reduction in the sorting size to the 100 nm and potentially below, enabling separation for biomedical analysis (e.g. exosomes)
- Tunability over a wide size range, the maximal sorting size of the device remains the same as for round shapes (micrometers, depending on the exact layout), while the minimal sorting size is less than 100 nm

Applications

- Exosome purification
- Sample purification for medical workflows
- Filtering of nanosized biological and chemical samples
- Water purification
- Quality control in pharma, cosmetics or chemistry

• Tunable size-based sorting